

PURETONE AUDIMED TINNITUS SYSTEMS


MM range of generic sound generators



Device	Type	Colour	Volume Control	Output Trimmer	Tone Trimmer	Vented	Battery Type	Fitting
MM1	Straight-stepped generic shell in-the-canal	White Blue Red Green	✓ ✓ ✓ ✓	✓		✓ ✓ ✓ ✓	10A	Unifit
MM2	Straight-contoured generic shell in-the-canal	White Blue Red Green	✓ ✓ ✓ ✓	✓		✓ ✓ ✓ ✓	10A	Unifit
MM4	Open ear generic in-the-helix	White Blue Red	✓ ✓ ✓	✓	✓ ✓		10A	Unifit
MM6 <small>Also Available - MM6D ((DIGITAL))</small>	Open ear generic in-the-helix	White Blue Red	✓ ✓ ✓	✓			10A	Left & Right Models
MM8 (C2008M)	Tapered generic shell in-the-canal	Blue	✓			✓	10A	Unifit
MM10 <small>Also Available - MM10D ((DIGITAL))</small>	Straight shell generic micro in-the-canal	White Blue Red Green	✓ ✓ ✓ ✓	✓	✓ ✓ ✓	✓ ✓ ✓ ✓	10A	Unifit

Output Levels

- White: 50dB
- Blue: 60dB
- Red: 80dB
- Green: 95dB



Battery Life

- White: 90 hours
- Blue: 60 hours
- Red: 60 hours
- Green: 60 hours

Puretone's policy is one of continuous development and improvement which may cause the data herein to vary without prior notice.

Puretone Ltd, 9-10 Henley Business Park, Trident Close, Medway City Estate, Rochester, Kent, UK.
T: +44 (0) 1634 719427 F: +44 (0) 1634 719450 E: info@puretone.net W: www.puretone.net

