

Series 1200 Technology



SoundScape technology lets you immerse yourself in the best available 360° listening experience, mimicking the organic sound processing of the human ear to produce sounds that are rich, clear and natural as they can be. Additionally, when you need it to, it automatically lets you focus on speech even in noisy surroundings.



2.4 GHz technology allows adjustments to be made wirelessly. Optional wireless accessories are also available for every listening experience.



The Feedback Manager Plus feature ensures feedback free hearing when you get close to someone or hold a phone close to your ear. Feedback Manager Plus targets and removes feedback, rapidly and precisely, without influencing the settings of the hearing device.



The Personal Noise Reduction feature automatically reduces noise depending on where you are. By a smaller degree when you're in the comfort of your own home, and by a much larger degree if you are in a noisy area, such as in a restaurant with family and friends.



The P Series state-of-the-art automatic directionality feature adjust to the sounds around you so you'll be able to focus on speech at all times.

Available from:



Series 1200 Hearing Systems

Advanced Digital Hearing Systems





Series 1200 Hearing Systems



Series 1200

Hearing solutions packed full with the latest technologies

Series 1200 hearing systems boast a wide range of features, offering a first class hearing experience. The range features Speech focused directionality while 2.4 GHz technology allows you to use the latest wireless accessories with your device.

The range includes 5 custom-made devices, as well as three behind-the-ear devices, with the P1200MW & P1200PW featuring optional open fittings and the P1200RW featuring receiver-in-the-ear technology.



 Wonder Ear™ 1200S ● ● Wireless ● 13	 Wonder Ear™ 1200P ● ● Wireless ● 13	
 CE2000™ 1200S ● ● ● Wireless ● 312	 CE2000™ 1200P ● ● ● Wireless ● 312	 CIC2010 1200P ● ● ● ● 10A
 P1200SW ● ● Wireless ● 13	 P1200MW ● ● Wireless ● 312	
 P1200PW ● ● Wireless ● 13	 P1200RW ● ● Wireless ● 312	