



# Maintenance & Hygiene

**Wax Removal Sprays**

**Conditioning Units**

**Ear Lotions**

**Earmould Cleaning Solutions**

**Brushes & Tools**



# Maintenance & Hygiene

## PerfectDry® Lux

UV-C disinfection and drying solution for hearing systems. Features forced-air fan drying system with 360 degree UV-C lamp. Powered via mains or USB. Supplied with UK mains adaptor.

NOW DRIES IN  
**30 MINS**



**ACC430H** PerfectDry® Lux  
**ACC430K** EU Power Supply  
**ACC430M** UK Power Supply

## PerfectDry®

Electric drying solution for hearing systems. Can hold 2 BTEs with earmoulds, several ITEs, glasses etc. Features secure drying cycle with timer and automatic switch-off. Effectively dries hearing systems without condensation. Supplied with UK mains adaptor.



**ACC430G** PerfectDry®  
**ACC430J** EU Power Supply

## Dry & Store®

Dry  
& Store

A complete conditioning system for all types of hearing instruments and in-the-ear devices. While you sleep, it removes moisture, dries earwax, kills germs and deodorises your hearing instrument. Dimensions: 180mm x 99mm x 91mm. Supplied with UK mains adaptor and 1 x desiccant Dry-Brik II.



**ACC357A** Dry & Store®  
**ACC357D** UK Mains Adaptor  
**ACC357H** Euro Mains Adaptor  
**ACC357G** 3 x Replacement Dry Brik II  
**ACC357G-3** 9 x Replacement Dry Brik II  
**ACC357G-10** 30 x Replacement Dry Brik II  
**ACC357C** Replacement Lamp

## Dry & Store® Zephyr

Dry  
& Store

Miniature automatic drying and conditioning system. Removes moisture from hearing instruments overnight. Dimensions: 150mm x 92mm x 72mm. Supplied with UK mains adaptor and 1 x desiccant Dry-Brik II.



**ACC357J** Dry & Store® Zephyr  
**ACC357G** 3 x Replacement Dry Brik II  
**ACC357G-3** 9 x Replacement Dry Brik II  
**ACC357G-10** 30 x Replacement Dry Brik II

## PerfectClean

The first system to clean, disinfect, and dry hearing aids with a press of a button

- Vortex cleaning action effectively removes earwax from hearing aid devices
- Helps to prevent irritation
- Maintains optimal performance of hearing aid devices
- Easy to use
- Fits all types (standard and custom-made)
- 3 actions over 90 minutes
  1. Cleaning: earwax removal solution
  2. Disinfecting: active anti-bacterial agent + UV-C system
  3. Drying: forced-air fan drying system
- Easy to refill
- Visual indication of the solution level



PerfectClean  
Refill



- ACC430P** PerfectClean
- ACC430Q** PerfectClean Refill
- ACC430R** PerfectClean EU Power Supply
- ACC430S** PerfectClean UK Power Supply

## Dry Spot Dehumidifier

Stylish container with microwaveable desiccant pillow. Designed to dry earwax and kill germs in hearing systems. Improves hearing system reliability and reduces the need for repairs.



- ACC333A** 1 x Dry Spot Dehumidifier
- ACC333A-10** 10 x Dry Spot Dehumidifier
- ACC333B** 1 x Spare Desiccant Pillow
- ACC333B-20** 20 x Spare Desiccant Pillow

## Suction Tool

Ideal for the removal of moisture and wax from hearing instruments. The housing separates to provide easy access for cleaning.



**ACC336** Suction Tool

## Moisture Control Kit

For drying hearing systems overnight. Contains: static shield bag, desiccant and 3-spot humidity card. Available as a kit or individually.



- ACC171A** 1 x Complete Kit
- ACC171A-50** 50 x Complete Kit
- ACC171B** 1 x Static Shield Bag
- ACC171C** 1 x Desiccant
- ACC171D** 1 x 3-Spot Humidity Card

## PowerVac Tips

Replacement vacuum tips for PowerVac system.



**ACC347B** Suction Tips

# Maintenance & Hygiene

## Earigator™ Cerumen Management System

Save time, retain your patients, and grow your practice with the Earigator™, the world's most advanced self-contained ear irrigation system. Full temperature and pressure control allows you to safely and comfortably clean out your patient's ear wax in under 30 seconds.

Designed by an Otologist, the Earigator™ combines the functions of an otoscope and irrigation into one, providing the most advanced means of cerumen control.

With the Earigator™, nurses and clinicians are fully capable of safely and efficiently removing even the most stubborn ear wax.

Self-contained temperature control constantly regulates the water to body temperature, avoiding any caloric or vertigo side effects.

Pressure controls ensure that the Earigator™ quickly removes even impacted cerumen, while never endangering the ear drum



Stainless Steel Cart



Suggested Buy



New Product

- ACC485A** Earigator
- ACC485B** Earigator Stainless Steel Cart
- ACC485C** 100 x Disposable Nozzles

## Ultraspec Single-Use Ear-Irrigation Noots Tank

New, improved noots-type ear tank facilitates easier, cleaner treatment of cerumen. The anti-splash tab significantly reduces back-wash during ear-irrigation procedures.

- Stable
- Leak resistant
- Space saving
- Eco-friendly

- ACC483-25** 25 Pack
- ACC483-50** 50 Pack
- ACC483-100** 100 Pack



New Product


# Maintenance & Hygiene

## Air Puffer

For removing moisture from tubing.  
Available in blue.



**Suggested Buy**



**ACC121** 1 x Air Puffer  
**ACC121-10** 10 x Air Puffer

## Acu-Life® Air Puffer

**ACU-LIFE®**

Helps keep earpieces and tubing clear from moisture, dirt and debris for clear, unobstructed sound. Great for earpieces, tubes and open fit hearing aids.



**ACC435C** Air Puffer

## Propulse® Ear Irrigator Pack

**PROPULSE®**  
IRRIGATOR™

Featuring a bold, modern design, the new Propulse® ear Irrigator sets new standards in ear irrigation whilst retaining the familiar and trusted features that you and your patients depend upon. The Propulse® ear Irrigator has been designed with both the patient and the clinician in mind. Ergonomically designed, the new handle design makes it easier to remove the QrX™ tips. Tried and tested key Propulse® features such as easy hands-free operation via the footswitch are maintained. A range of products & accessories that complement the Propulse® ear Irrigator are also available.

- Trusted Propulse® quality
- Hands-free operation by use of footswitch
- Variable pressure control for comfort and safety
- Ergonomically designed handle
- Easy tip removal & replacement
- Utilises Propulse® QrX™ Technology



**ACC458A** Propulse® Ear Irrigator Pack

## Propulse® QRX Tips

**PROPULSE®**

Propulse® Single Use QrX™ Tips are designed exclusively for use with the Propulse® G5 Ear Irrigator. The unique locking mechanism provides an easy fit on to the quick release handle whilst securing the tip throughout the irrigation procedure. Supplied in pack of 100.



**ACC458B** 100 x QRX™ Tips - White

## Propulse® ProTect™

**PROPULSE®**

This waterproof shoulder cape helps protect the patient's clothes during the ear irrigation procedure. Supplied in pack of 36.



**ACC458C** ProTect™

# Maintenance & Hygiene

## Propulse® Mushroom Valve & Washer

The Mushroom Valve is specifically designed to prevent water flowing out of the reservoir whilst filling.



ACC458D Mushroom Valve & Washer

## Propulse® Carry Case

- Soft, padded carry case designed to carry the Propulse®
- Zipped flap allows quick and easy access to the Propulse®
- Prevents damage to the Propulse® during transit
- Suitable for domiciliary visits
- Convenient shoulder strap



ACC458G Carry Case

## Propulse® Cleaning Tablets

Propulse® Cleaning Tablets are specially formulated for use with Propulse® Ear Irrigators.

They ensure that the irrigator is cleaned effectively without damaging the pump and working components of the unit as some other cleaning agents can. Supplied in pack containing 200 tablets.



ACC459A Cleaning Tablets

## Propulse® ProScoop

The double ended design combines a scoop for manual ear wax removal at one end, and a serrated edge at the other to allow easy attachment of cotton wool (for mopping out the ear following ear irrigation). Supplied in pack of 50.



ACC460A ProScoop

## Propulse® Head Lamp

Illumination of the ear canal is essential when performing ear irrigation.

- Flexible goose neck design provides flexibility of movement
- Comfortable and light weight
- Bright white light provides excellent illumination of the external auditory meatus
- Elasticated headband fits all users
- Simple push button on/off switch



ACC461A Propulse® Head Lamp

ACC461B Propulse® Head Lamp Bulb

## Propulse® Noots Ear Tank

The Propulse® reusable noots ear tank is ideal for ear irrigation and other cerumen removal procedures.



ACC462A Noots Ear Tank

## Acu-Life® Ear Irrigator **ACU-LIFE®**

The Acu-Life® ear irrigator cleaning system cleans ears and helps to prevent ear wax buildup. Tri-stream soft tip and pressure control bottle to ensure safe and effective cleaning. Bottle contains 118ml.



**ACC435D** Ear Irrigator

## Acu-Life® Wax Removal Syringe **ACU-LIFE®**

Ear wax removal 20ml syringe for the use with water or ear wax drops. Safely clean ears with the soft tri-stream tip.



**ACC435B** Wax Removal Syringe

## Acu-Life® Oto-Scoop **ACU-LIFE®**

Oto-scoop with flexing tip is the safe and simple way to remove unwanted ear wax. Simply insert flexing tip into the ear and sweep the ear wax out. Safety guard protects ear drum. For adult use only.



**ACC435F** Oto-Scoop

## Wax Remover Sprays

For effective removal of earwax from hearing instruments and in-ear monitors. Available in 50ml and 30ml sizes (please note: 50ml size features compressed gas - special delivery only).



**Special Delivery**

50ml  
Compressed Gas  
Hazardous Shipping

- ACC309A** PureClean (50ml)
- ACC309A-3** 3 x PureClean (50ml)
- ACC309A-24** 24 x PureClean (50ml)



**Special Delivery**

50ml  
Compressed Gas  
Hazardous Shipping

- ACC309B** IEM Spray (50ml)
- ACC448G** PureClean (30ml)
- ACC448G-10** 10 x PureClean (30ml)



**Suggested Buy**

30ml  
Pump Action

# Maintenance & Hygiene

## PureClean Cleaning Wipes

Easy to use moist, hygienic cleaning wipes. Ideal for a quick cleaning or use while travelling. Supplied individually wrapped in boxes of 30 wipes.



**Suggested Buy**

**ACC448A** Single Pack (30 Wipes)  
**ACC448A-10** 10 Packs (300 Wipes)

## PureClean Comfort Gel

Earmould sealant/comfort gel for easier insertion and removal of in-ear hearing systems, earpieces and earmoulds. Supplied in 5ml tube. Supplied individually or in a pack of 10.



**Suggested Buy**

**ACC448B** Single Pack  
**ACC448B-10** 10 Pack

## PureClean Cleaning Tablets & Cup

Used regularly they can easily remove earwax accumulated through normal wear. Tablets supplied in pack of 20 or 200. Cleaning cups supplied individually or as packs of 10.



**ACC448E** 1 x Cleaning Cup  
**ACC448E-10** 10 x Cleaning Cups  
**ACC448F** 20 x Cleaning Tablets  
**ACC448F-10** 200 x Cleaning Tablets

## PureClean Drying Capsules & Cup

Ideal for everyday use, they absorb any moisture present in hearing systems. Drying Capsules supplied in pack of 4 or 40. Drying cups supplied individually or as packs of 10.



**ACC448C** 1 x Drying Cup  
**ACC448C-10** 10 x Drying Cups  
**ACC448D** 4 x Drying Capsules  
**ACC448D-10** 40 x Drying Capsules

## PureClean Hearing Aid Starter Set

Contains 6 PureClean cleaning wipes, a PureClean cleaning spray, a PureClean drying capsule and a PureClean drying cup. Includes zip bag, brochure and PureClean hearing aids sleeve.



**ACC447A** Hearing Aid Starter Set

## PureClean Noise Protection Starter Set

Contains 6 PureClean cleaning wipes, a PureClean comfort gel and an air puffer. Includes zip bag, brochure and PureClean noise protection sleeve.



**ACC447B** Noise Protection Starter Set



## Maintenance Pack

All-in-one hearing system maintenance kit. Supplied with zippa bag, dehumidifier bag, stetoclip, tubing extractor, tubing threader, blue air puffer, PureClean hygienic wipes and 5 x pre-bent tubes.



**ACC302** Maintenance Pack

## Medi Wipes


Alcohol based cleaning wipes for hearing systems and earmoulds. Available in a box of 100, 1000 and individually.



- ACC404-1** Individual
- ACC404** Pack of 100
- ACC404-10** Pack of 1000

## Cleaning Tool Kits

Convenient pocket sized multi-tools for hearing systems. Features magnet, battery tray opener, vent cleaner, wax loop, cleaning brush, as well as 2 battery storage compartments. Available with or without digital battery tester.



- ACC398** 1 x Without Battery Tester
- ACC398-10** 10 x Without Battery Tester
- ACC421** 1 x With Battery Tester
- ACC421-5** 5 x With Battery Tester

## Wonder Wipes


Anti-bacterial cleaning wipes for hands, tools and surfaces. Resistant to MRSA, Salmonella, Listeria & E. Coli. Supplied in tub of 100 wipes.



- ACC431** Single Tub
- ACC431-3** 3 x Tubs

## Multi-Bristle Vent Cleaner

Vent cleaning tool with a selection of multi-sized bristles.



- ACC392** Multi-Bristle Vent Cleaner

# Maintenance & Hygiene

## AudiSpray®

100% natural solution for ear hygiene, Audispray® gently cleanses the ears by dispersing the earwax. Audispray® for age 6 and above, and is supplied in a 50ml bottle (about 6 months of use).



**ACC424A** Adult  
**ACC424A-PACK** 9 x Adult  
**ACC424C** Display Box

## AudiClean®

AudiClean® spray is the safer alternative to cotton buds. Made from purified seawater, AudiClean® spray gently cleanses the ears by dispersing the earwax.

Used on a regular basis AudiClean® spray will help prevent the formation of earwax. Contents 115ml.



## Eargene®

Soothing and refreshing ear lotion. Relieves itching ears, leaving the ear cool and refreshed. Supplied in a 15ml bottle, available individually or in a display pack of 24.



**Suggested Buy**

**ACC286** 1 x 15ml  
**ACC286-24** 24 x 15ml

**ACC375E** Single  
**ACC375E-10** Pack of 10

## Earol® Olive Oil Spray

Earol® Spray has a simple yet unique delivery system capable of dispensing a discrete unit dose of olive oil into the outer ear, softening compacted earwax. Earol® can be used preventively once a week to aid ear hygiene and maintain healthy ears. Supplied in 10ml spray.



**Suggested Buy**

**ACC449** 1 x Earol  
**ACC449-10** 10 x Earol

# Maintenance & Hygiene

## Brushes & Tools

Wax Loop (ACC114)



Power Magnet & Wax Loop (ACC117J)



Brush & Wax Loop (ACC117C)



Brush (ACC116)



Wax Loop & Magnet (ACC115)



CENS ProFlex Brush & Magnet (ACC117G)



Brush & Magnet (ACC117A)



Brush & Child Lock Tool (ACC116L)



Wonder Tool Battery Remover (ACC213)



Brush & Vent Cleaner (ACC117D)



Brush, Magnet & Wax Loop (ACC117B)



Anticer Brush (AN3)



**ACC114** Wax Loop

**ACC115** Wax Loop & Magnet

**ACC116** Brush

**ACC116L** Brush & Child Lock Tool

**ACC117A** Brush & Magnet

**ACC117B** Brush, Magnet & Wax Loop

**ACC117C** Brush & Wax Loop

**ACC117D** Brush & Vent Cleaner

**ACC117G** CENS® ProFlex Brush & Magnet

**ACC117J** Power Magnet & Wax Loop

**ACC213** Wonder Tool

**AN3** Anticer Brush

**ACC114-150** 150 x Wax Loops

**ACC115-100** 100 x Wax Loop & Magnet

**ACC116-100** 100 x Brushes

**ACC116L-100** 100 x Brush & Child Lock Tools

**ACC117A-100** 100 x Brush & Magnets

**ACC117B-100** 100 x Brush, Magnet & Wax Loops

**ACC117C-100** 100 x Brush & Wax Loops

**ACC117D-100** 100 x Brush & Vent Cleaners

**ACC117G-100** 100 x CENS® ProFlex Brush & Magnets

**ACC117J-50** 50 x Power Magnet & Wax Loops

**ACC213-10** 10 x Wonder Tools

**AN3-50** 50 x Anticer Brushes

# Maintenance & Hygiene

## PureVac

A revolutionary hearing aid vacuum system featuring a stylish aluminium enclosure, 25 bar vacuum pressure, and an all-in-one, compact design. This sleek, modern system is perfect for any hearing care professional.

Compact, modern and simple to use device for the cleaning of hearing instruments. Individual vacuum and pressure styluses allow for hygienic concentrated cleaning of a device. A timed drying chamber utilising both vacuum and conducted heat provides all over cleansing.

- 25 bar vacuum
- 50°C typical chamber heat
- 7 minute run time
- Compact design
- Easily replaceable tube
- Easily replaceable stylus
- Pressure display



## Pure Vac



Suggested Buy

ACC473 PureVac

## Jodi-XL

Complete clinical vacuum and desiccator system in a small, yet powerful pump. Two 2000 hour-rated vacuum pumps power this unit; one pump draws moisture and debris out through a needle vacuum and the other powers the heated desiccator. Includes a selection of different-sized needles for vents and sound bores. Powered by a 15V DC transformer, this unit can be used anywhere in the world.



ACC411A Jodi-XL

ACC411B Filters - Pack of 3

ACC411C Filters - Pack of 6