



Tinnitus Therapy

Sleep Sound Therapy
Tinnitus Relief Devices
Pillow Speakers



Tinnitus Therapy

Sound Oasis® (V3)


The Sound Oasis® bedside desktop sound generator & alarm clock features 12 natural sounds to soothe and relieve the symptoms of tinnitus at night:

- Ocean Surf
- Rain
- Stream
- White Noise
- Woodlands
- Thunderstorm
- Energy Chimes
- Tranquillity
- Alpha Clouds
- Sleep Surf
- Delta Voyage
- Dream Echoes

Features a built-in timer with settings for 30, 60 or 90 minutes, a built-in alarm clock with snooze button, as well as an LCD display. The sound generator can be powered either from the UK mains adaptor (supplied) or 4 x AA batteries (available separately). Ideal for use with the Sound Pillow (ACC386). Additional sound cards also available (see below).



- ACC385E** Sound Oasis® (V3)
- ACC385E-2** 2 x Sound Oasis® (V3)
- ACC385FV2** Sound Oasis® (V2) UK Power Supply
- ACC385FV3** Sound Oasis® (V3) UK Power Supply
- ACC385C** EU Mains Adaptor



Sound Oasis® Sound Cards

Sound Oasis® Sound Cards allow Sound Oasis® owners to expand their collection of available sounds. By simply inserting the sound card into slot located on the base of the Sound Oasis®, users can supplement their unit's built-in sound offerings with 12 additional digitally recorded sounds (compatible with Sound Oasis® V3 only). Supplied as either a single card or a pack of 5.

Ear Therapy - 12 Tracks: Ocean surf, Rain, Stream, White noise 4k, White noise 6k, Full spectrum, Thera1 1k-10k, Thera2 1k-10k, Thera3 1k-10k, Thera 2k-8k, Thera 2.5k-8k, Thera 2.5k-5k, Thera .9k-3.2k

Spa Retreat - 12 Tracks: Hawaiian surf, Gentle rain, Waterfall, Island stream, Asian garden, Summer night, Unwind, Reflection, Serenity, Easy days, Late night, Stillness

Nature Journey - 12 Tracks: Thunderstorm, Crashing waves, Forest rain, Waterfall, Gentle surf, Ocean rain, Wind, Secluded lake, Caribbean surf, Jungle stream, Summer night, Amazon lagoon

Sleep Sounds For Baby - 12 Tracks: White noise, Ocean surf, Gentle rain, Summer night, Dolphins, Prenatal, Heartbeat, Womb, Lullaby, Car ride, Lullaby 2, Woodlands

- ACC385-C8** Ear Therapy
- ACC385-C8-5** 5 x Ear Therapy
- ACC385-C9** Spa Retreat
- ACC385-C9-5** 5 x Spa Retreat
- ACC385-C10** Nature Journey
- ACC385-C10-5** 5 x Nature Journey
- ACC385-C11** Sleep Sounds For Baby
- ACC385-C11-5** 5 x Sleep Sounds For Baby



Sound Oasis® Mini Sound Machine



New Product

The World's Smallest Sound Machine® will help you relax, fall asleep easier, get higher quality sleep, manage tinnitus, and even improve concentration during the day. Our unsurpassed collection of therapeutic sounds are selected to provide you with the best possible solution for your needs. 15 hour battery life.



ACC385Q Mini Sound Machine

Sound Oasis® Mini White Noise Machine



New Product

Enjoy soothing white noise wherever you go! The latest innovation from Sound Oasis® gives you white noise in its most portable form yet. Up to 20 hours of run time per charge.



ACC385P Mini White Noise Machine

Sound Oasis® Deluxe

The Sound Oasis® Deluxe bedside desktop sound generator & alarm clock features a massive 145 authentic nature sounds to soothe and relieve the symptoms of tinnitus at night.

Features professional sound quality, thanks to its 3 speaker (including subwoofer), and also features an AM/FM radio, dual 12/24 hour clock with 3 alarms and a digital calendar. Includes 4 position timer and memory button, as well as a headphone jack and a CD/MP3 input jack. Powered by AC adaptor (included).



ACC385G Sound Oasis® Deluxe

Sound Oasis® Bluetooth™

Sleep better with any one of 10 built-in sounds or stream via Bluetooth™ from your mobile device or computer. Perfect for improving sleep, blocking noise, and managing tinnitus. Features high quality speaker for clear sound quality with powerful bass. Powered by built-in rechargeable battery for wireless freedom or from USB. 30'/10m Bluetooth™ wireless range.



ACC385L Sound Oasis® Bluetooth™

Sound Oasis® S-550-06

Natural sounds are a proven, healthy sleep aid. By bathing your room with nature sounds, the S-550-06 sleep machine provides effective, natural, prescription free sleep therapy.

Ideal for help with sleep or managing tinnitus. For extra sound variety, add more sounds with the patented sound card system.



ACC385M Sound Oasis® S-550-06

Sound Oasis® Travel

The Sound Oasis® Travel bedside desktop sound generator & alarm clock features 19 authentic nature sounds to soothe and relieve the symptoms of tinnitus at night.

Features a backlit, dual 12/24 hour alarm clock with a digital calendar. Includes 4 position timer and memory button, as well as a headphone jack. Powered by AC adaptor (included). Travel pouch and time zone chart included.



ACC385H Sound Oasis® Travel

ACC385H-2 2 x Sound Oasis® Travel

Tinnitus Therapy

Sound Oasis® Pillow Speakers

Sleep while listening to your favourite sounds!
Turn any pillow into a Sleep Therapy Pillow with these ultra-thin stereo speakers. Simply place under your existing pillow and enjoy your favourite music or sounds for optimal relaxation and comfort. With the option of built-in amplification and two high fidelity, ultra-thin stereo speakers, the amplified pillow speaker is perfect for surrounding one's head with music and sound for enhanced therapy or privacy. Ideal for help with sleep or managing tinnitus.



- PS2** Stereo Pillow Speaker
- PS3** Amplified Stereo Pillow Speakers

Sound Pillow®

The Sound Pillow® looks and feels just like a normal pillow, but features built-in stereo speakers so that the sound from a bedside sound generator can be discreet and intimate. The Sound Pillow® is fitted with a cable and 3.5mm jack plug for plugging into the Sound Oasis® sound generator, a personal stereo or clock radio. Available with or without volume control. (Not compatible with the Naturcare® Sounds Relaxation Ball).
Sound Pillow®: 59cm x 39cm x 16cm. Travel Sound Pillow®: 39cm x 24cm x 10cm.



Please note: This product is non-returnable.

- ACC386B** Sound Pillow®
- ACC386B-3** 3 x Sound Pillow®
- ACC386C** Travel Sound Pillow®

Naturcare® Sounds Relaxation Ball

Aids tinnitus sufferers in getting a restful night's sleep by playing digital recordings of real, natural sounds. Choice of 7 digitally recorded sounds: Heartbeat, Sea Waves, Brook, Birds, Woods, White Noise, Rain.

Features a volume control, a 30 minute timer, as well as no tapes or moving parts to wear out. 4 x AA batteries included, with UK and EU power adaptors available separately. Bundle pack with UK adaptor also available.

- ACC323** Naturcare® Sounds Relaxation Ball
- ACC323A** Naturecare® Bundle Pack
- ACC323B** UK Mains Adaptor
- ACC323B-10** 10 x UK Mains Adaptors
- ACC323C** EU Mains Adaptor



Suggested Buy

Tinnitus Therapy

GTS-2000 DELUXE Glo to Sleep Therapy Mask

DELUXE Glo to Sleep is a safe, natural and effective way to help you get the sleep you need. Ancient wisdom joins with modern technology creating a revolutionary new way for you to clear your mind, quickly relax and get to sleep naturally. Sound Oasis® DELUXE Glo to Sleep will help you focus your mind on the process of relaxation and sleep - naturally and effectively. Features adjustable brightness with 3 auto dim options: ramp down in 10, 20 or 30 minutes.



ACC385J GTS-2000 DELUXE

Tinnitus Relief Devices

PURETONE
AUDIMED
TINNITUS SYSTEMS



MM1 Generic Sound Generator

Generic straight stepped in-the-canal white noise generator for tinnitus. Features volume control, vent and type 10 battery (50dB model also features output trimmer).

- MM1-WHITE** MM1 - 50dB Output
- MM1-BLUE** MM1 - 60dB Output
- MM1-RED** MM1 - 80dB Output
- MM1-GREEN** MM1 - 95dB Output

MM6 Generic Sound Generator

Generic in-the-helix white noise generator for tinnitus. Features volume control and type 10 battery (50dB model also features output trimmer). Available in left (L) and right (R) ear models.

- MM6-WHITE-L** MM6 - 50dB Output - L Ear
- MM6-WHITE-R** MM6 - 50dB Output - R Ear
- MM6-BLUE-L** MM6 - 60dB Output - L Ear
- MM6-BLUE-R** MM6 - 60dB Output - R Ear
- MM6-RED-L** MM6 - 80dB Output - L Ear
- MM6-RED-R** MM6 - 80dB Output - R Ear

MM2 Generic Sound Generator

Generic straight contoured in-the-canal white noise generator for tinnitus. Features volume control, vent and type 10 battery (50dB model also features output trimmer).

- MM2-WHITE** MM2 - 50dB Output
- MM2-BLUE** MM2 - 60dB Output
- MM2-RED** MM2 - 80dB Output
- MM2-GREEN** MM2 - 95dB Output

MM8+ Generic Sound Generator

Generic tapered shell, in-the-canal white noise generator for tinnitus. Features volume control and vent.

- MM8+** MM8+ - 60dB Output

MM10 Generic Sound Generator

Generic ultra-miniature straight contoured in-the-canal white noise generator for tinnitus. Features volume control, vent and type 10 battery (50dB model also features output trimmer).

- MM10-WHITE** MM10 - 50dB Output
- MM10-BLUE** MM10 - 60dB Output
- MM10-RED** MM10 - 80dB Output
- MM10-GREEN** MM10 - 95dB Output

Tinnitus Therapy

MM10-DAC TRT

Digitally programmable (via Purefit 5) generic ultra-miniature straight contoured in-the-canal white noise generator for tinnitus. Features volume control, vent and type 10 battery.

PURETONE
AUDIMED
TINNITUS SYSTEMS



MM10-DAC TRT Beige

MM6-DAC TRT

Digitally programmable (via Purefit 5) generic in-the-helix white noise generator for tinnitus. Features volume control and type 10 battery. Available in left and right ear models.

PURETONE
AUDIMED
TINNITUS SYSTEMS



MM6-DAC TRT-L Left Ear
MM6-DAC TRT-R Right Ear

MM6 Tubing Assembly

Replacement tubing assembly for mm6 and mm6-DAC TRT models. Supplied individually or in packs of 20.



MM6TUBE-ASSY Single
MM6TUBE-ASSY-20 20 Pack

22MA+ BTE

Miniature white noise generating BTE, featuring volume control, on/off switch and a type 312 battery. Available in beige or brown. Also available as open fit.

PURETONE
AUDIMED
TINNITUS SYSTEMS



Suggested Buy

22MA+BEI Beige
22MA+OPEN-BEI Beige - Open

Anti-Feedback Sleeve

Ideal for CIC hearing systems, for eliminating feedback in the ear. Made of thin, soft silicone which can be cut down to size as needed. Also used with themm range of tinnitus devices. Supplied individually or in packs of 50.



ACC168 Single
ACC168-50 50 Pack

PT-40 DAC TRT

Miniature white noise generating BTE, featuring volume control, on/off switch and a type 312 battery. Digitally programmable via Purefit 5. Available in beige.

PURETONE
AUDIMED
TINNITUS SYSTEMS



PT-40-DAC TRT-BEI Beige