



Tubing, Tips, Earmoulds & Hooks

Earmould Tubing
Hearing System Tubing
Tubing Tools
General Purpose Ear Tips
Generic Earmoulds
BTE Earhooks

Tubing Size Guide

Size	Inside Diameter	Outside Diameter	Wall Thickness
Standard	2.160mm	3.170mm	0.505mm
Medium	1.930mm	3.100mm	0.580mm
Medium Moisture Free	1.930mm	3.100mm	0.580mm
Thickwall	1.930mm	3.600mm	0.830mm
Thickwall Moisture Free	1.930mm	3.300mm	0.680mm



Tubing, Tips, Earmoulds & Hooks

Pre-Bent Tubing


Our most popular tubing, these pre-bent lengths of tubing help save time and effort when fitting BTE hearing systems. Available in 3 sizes and supplied in packs of 25, 250 & 500. Moisture free variation available in 2 sizes and supplied in packs of 5 & 50.



ACC165A	25 x Standard
ACC165A-10	250 x Standard
ACC165A-20	500 x Standard
ACC165E	25 x Medium
ACC165E-10	250 x Medium
ACC165E-20	500 x Medium
ACC165B	25 x Thickwall
ACC165B-10	250 x Thickwall
ACC165B-20	500 x Thickwall
ACC165F	5 x Moisture Free - Medium
ACC165F-50	50 x Moisture Free - Medium
ACC165G	5 x Moisture Free - Thickwall
ACC165G-50	50 x Moisture Free - Thickwall

Single Bend Tubing

Single bend tubing available in medium and thickwall sizes. Supplied in packs of 25 tubes.



ACC165E-SB	25 x Medium
ACC165B-SB	25 x Thickwall

DisappEar Pre-Bent Tubing

The smooth, matte finish of DisappEar tubing makes it very discreet. Available in standard, medium and thickwall sizes. Supplied in packs of 10 tubes.



ACC165H	10 x Standard
ACC165J	10 x Medium
ACC165K	10 x Thickwall

Pre-Bent Tubing With Tube Lock

Pre-bent tubing that uses a brass retention ring for soft earpieces.



ACC160K	10 x Tubing With Tube Lock
----------------	----------------------------

Libby Horn Tubing

These pre-bent, internally stepped, tapered lengths of tubing help preserve and extend high frequency response. Available in standard and moisture free variations, in 3mm and 4mm sizes. Supplied individually or in packs of 10.



ACC221A	1 x 3mm
ACC221A-10	10 x 3mm
ACC221B	1 x 4mm
ACC221B-10	10 x 4mm
ACC221E	1 x Moisture Free - 3mm
ACC221E-10	10 x Moisture Free - 3mm
ACC221F	1 x Moisture Free - 4mm
ACC221F-10	10 x Moisture Free - 4mm

Tubing, Tips, Earmoulds & Hooks

Link Tubing

Available in 4 sizes. Supplied per metre or on 305m/1000ft roll (Super Thickwall only available in a roll).



ACC163	Standard
ACC170	Medium
ACC164	Thickwall
ACC163ROLL	Standard Roll
ACC170ROLL	Medium Roll
ACC164ROLL	Thickwall Roll
ACC169ROLL	Super Thickwall Roll

Receiver Tubing

Outside diameter 2.1mm, inside diameter 1.5mm. Available in 3 colours and supplied per metre. Pre-bent/pre-cut variation supplied individually.



ACC222D	Pink
ACC222A	Clear
ACC222Q	Pre-Bent/Pre-Cut - Black
ACC222U	Pre-Bent/Pre-Cut - Clear

90 Degree Bends

Outside diameter 3mm. Available in clear or pink, in 2 sizes. Supplied individually or in packs of 100.



ACC160A	1 x Clear - Short Shank
ACC160A-100	100 x Clear - Short Shanks
ACC160B	1 x Pink - Short Shank
ACC160B-100	100 x Pink - Short Shanks
ACC160C	1 x Clear - Long Shank
ACC160C-100	100 x Clear - Long Shanks
ACC160D	1 x Pink - Long Shank
ACC160D-100	100 x Pink - Long Shanks

Semi-Rigid Tubing

Outside diameter 2mm. Supplied per metre.



ACC165D	Semi-Rigid Tubing
----------------	-------------------

Curved Bend

Hard plastic curved bend with outside diameter of 3mm. Supplied individually or in packs of 100.



ACC160E	Single
ACC160E-100	100 Pack

Tubing, Tips, Earmoulds & Hooks

Microphone Tubing


For ITE and ITC hearing systems. Outside diameter 1.6mm, inside diameter 1mm. Supplied per metre.



ACC222B Microphone Tubing

Transducer Tubing

Pink transducer tubing. Outside diameter 1mm. Supplied per metre.



ACC222C Transducer Tubing

Receiver Adaptor

Available in clear and pink.
Supplied individually or in packs of 100.



- ACC145A** 1 x Clear
- ACC145A-100** 100 x Clear
- ACC145B** 1 x Pink
- ACC145B-100** 100 x Pink

Ear Tip Adaptors

For connecting ear tips to body aid receivers.
Supplied individually or in packs of 20. Available in long and short.



- ACC402A** 1 x Short
- ACC402A-20** 20 x Short
- ACC402** 1 x Long
- ACC402-20** 20 x Long

Earmould Adaptor

Available in clear and pink.
Supplied individually or in packs of 100.



- ACC144A** 1 x Clear
- ACC144A-100** 100 x Clear
- ACC144B** 1 x Pink
- ACC144B-100** 100 x Pink

Tubing Terminators

For use with spectacle devices. Available with tubing (supplied in packs of 25) or without tubing (supplied individually).



- ACC221D** 25 x Tubed
- ACC221C** 1 x Non-Tubed

3-Way Tubing Adaptor

Supplied individually.



ACC145C Tubing Adaptor

Tubing, Tips, Earmoulds & Hooks

Puretone Easy Thin Tube Box

Puretone open fitting box, suitable for use with: P Series, PT-40, PT-70 and 22MA+. Contains various sizes of domes and tubes, universal tool, cleaning wire and a measuring tool.



Suggested Buy



PT-OP-FIT-KIT Easy Thin Tube Box

Retubing Kit

Kit contains: one pair of pliers, one tubing expander, a surgical blade, rounded end scissors, small screw driver and pointed tweezers. All presented in a high quality leatherette case.



ACC454 Retubing Kit

EzEFIT Earpieces



These anatomically shaped generic NanoComm earpieces can be used as an alternative to a custom earmould. The shaped barbed adaptor follows the bend of the ear canal and allows the user to fit either large or small dome soft silicone ear tips, providing a secure seal in the ear for excellent attenuation.

This also still allows for sound to be transferred into the ear from the headset for clear communications. Available individually or as a kit containing left and right earpieces in both sizes.



- NC-GT-ML Left - Medium
- NC-GT-LL Left - Large
- NC-GT-MR Right - Medium
- NC-GT-LR Right - Large
- NC-GT-KIT Earpiece Kit

Tubing Expander

Stainless steel tubing expander.



Suggested Buy



- ACC162B Single
- ACC162B-5 5 Pack

Blunt End Scissors

Ideal for trimming tubing quickly and safely.



Suggested Buy



- ACC307 Single
- ACC307-10 10 Pack

Tubing, Tips, Earmoulds & Hooks

Tubing Tools

These handy threader, extractor, and combination threader/extractor tools are ideal for quickly threading and removing tubing from earmoulds. Designed and manufactured by Puretone. Can be personalised (laser marked/inked). Please contact us for details.

Tubing Threader



Tubing Extractor



Tubing Extractor/ Threader



Suggested Buy

- ACC161** 1 x Tubing Threader
- ACC162C** 1 x Tubing Extractor
- ACC162D** 1 x Tubing Extractor & Threader

- ACC161-5** 5 x Tubing Threaders
- ACC162C-5** 5 x Tubing Extractors
- ACC162D-5** 5 x Tubing Extractor & Threaders

Comply Canal Tube Adaptors

COMPLY.

Supplied in packs of 15 or 45.



- ACC262B** 15 x Tubes
- ACC262B-3** 45 x Tubes

Comply Canal Tip Set

COMPLY.

This professional fitting system comprises of sound tube adaptors and slow recovery foam ear tips that can be squeezed so they expand in the ear to fit the natural contours of the ear canal. Includes 20 foams ear tips and 10 tube adaptors. Supplied in a presentation case.



Comply Tips

COMPLY.

Foam ear tips. Available in 2 sizes. Supplied in pairs.



- ACC262C** Small
- ACC262D** Large

- ACC262** Single
- ACC262-2** 2 Pack

Comply Canal Tips Refill Kits

COMPLY

Refill kit consisting of 12 slow recovery foam tips. Available in standard, slim and short with 0, 1, 2, 3 or 4 side vents. Supplied individually or in packs of 3.



ACC262A-10	1 x Standard - No Vent
ACC262A-10-3	3 x Standard - No Vent
ACC262A-11	1 x Standard - 1 Vent
ACC262A-11-3	3 x Standard - 1 Vent
ACC262A-12	1 x Standard - 2 Vents
ACC262A-12-3	3 x Standard - 2 Vents
ACC262A-13	1 x Standard - 3 Vents
ACC262A-13-3	3 x Standard - 3 Vents
ACC262A-14	1 x Standard - 4 Vents
ACC262A-14-3	3 x Standard - 4 Vents
ACC262A-20	1 x Slim - No Vent
ACC262A-20-3	3 x Slim - No Vent
ACC262A-21	1 x Slim - 1 Vent
ACC262A-21-3	3 x Slim - 1 Vent
ACC262A-22	1 x Slim - 2 Vents
ACC262A-22-3	3 x Slim - 2 Vents
ACC262A-23	1 x Slim - 3 Vents
ACC262A-23-3	3 x Slim - 3 Vents
ACC262A-24	1 x Slim - 4 Vents
ACC262A-24-3	3 x Slim - 4 Vents
ACC262A-30	1 x Short - No Vent
ACC262A-30-3	3 x Short - No Vent
ACC262A-31	1 x Short - 1 Vent
ACC262A-31-3	3 x Short - 1 Vent
ACC262A-32	1 x Short - 2 Vents
ACC262A-32-3	3 x Short - 2 Vents
ACC262A-33	1 x Short - 3 Vents
ACC262A-33-3	3 x Short - 3 Vents
ACC262A-34	1 x Short - 4 Vents
ACC262A-34-3	3 x Short - 4 Vents

Soft Shaped Ear Tips

Soft shaped clear ear tips. Available in 3 sizes, with or without tubes. Supplied with tubes individually or as packs of 10 or 100. Also available individually without tubes.



ACC102A	1 x Left - Small
ACC102A-10	10 x Left - Small
ACC102A-100	100 x Left - Small
ACC102D	1 x Right - Small
ACC102D-10	10 x Right - Small
ACC102D-100	100 x Right - Small
ACC102B	1 x Left - Medium
ACC102B-10	10 x Left - Medium
ACC102B-100	100 x Left - Medium
ACC102E	1 x Right - Medium
ACC102E-10	10 x Right - Medium
ACC102E-100	100 x Right - Medium
ACC102C	1 x Left - Large
ACC102C-10	10 x Left - Large
ACC102C-100	100 x Left - Large
ACC102F	1 x Right - Large
ACC102F-10	10 x Right - Large
ACC102F-100	100 x Right - Large
ACC102G	No Tube - Small
ACC102H	No Tube - Medium
ACC102J	No Tube - Large

Tubing, Tips, Earmoulds & Hooks

Anti-Feedback Soft Ear Tips

Available in multiple sizes in clear or pink, with or without tubes. Pink with tubes supplied either individually, in packs of 10 or 100. Clear with tubes supplied either individually or in packs of 100. Pink/clear without tubes supplied individually.

Pink With Tubes

ACC101A	1 x Pink - With Tube - Very Small
ACC101B	1 x Pink - With Tube - Small
ACC101C	1 x Pink - With Tube - Medium
ACC101D	1 x Pink - With Tube - Large
ACC101E	1 x Pink - With Tube - Very Large
ACC101F	1 x Pink - With Tube - Extra Large

ACC101A-10	10 x Pink - With Tube - Very Small
ACC101B-10	10 x Pink - With Tube - Small
ACC101C-10	10 x Pink - With Tube - Medium
ACC101D-10	10 x Pink - With Tube - Large
ACC101E-10	10 x Pink - With Tube - Very Large
ACC101F-10	10 x Pink - With Tube - Extra Large

ACC101A-100	100 x Pink - With Tube - Very Small
ACC101B-100	100 x Pink - With Tube - Small
ACC101C-100	100 x Pink - With Tube - Medium
ACC101D-100	100 x Pink - With Tube - Large
ACC101E-100	100 x Pink - With Tube - Very Large
ACC101F-100	100 x Pink - With Tube - Extra Large

Clear With Tubes


ACC101G	1 x Clear - With Tube - Very Small
ACC101H	1 x Clear - With Tube - Small
ACC101J	1 x Clear - With Tube - Medium
ACC101K	1 x Clear - With Tube - Large
ACC101L	1 x Clear - With Tube - Very Large
ACC101M	1 x Clear - With Tube - Extra Large

ACC101G-100	100 x Clear - With Tube - Very Small
ACC101H-100	100 x Clear - With Tube - Small
ACC101J-100	100 x Clear - With Tube - Medium
ACC101K-100	100 x Clear - With Tube - Large
ACC101L-100	100 x Clear - With Tube - Very Large
ACC101M-100	100 x Clear - With Tube - Extra Large

Pink/Clear Without Tubes

ACC101Q	Pink - No Tube - Very Small
ACC101N	Pink - No Tube - Small
ACC101P	Pink - No Tube - Large
ACC101R	Pink - No Tube - Very Large

ACC101S	Clear - No Tube - Very Small
ACC101T	Clear - No Tube - Small
ACC101U	Clear - No Tube - Medium
ACC101V	Clear - No Tube - Large
ACC101W	Clear - No Tube - Very Large



Tubing, Tips, Earmoulds & Hooks

Paediatric Ear Tips

3.5mm ear tips.
Available in red or blue.
Supplied individually or
in packs of 50.



ACC340A	1 x Red
ACC340A-50	50 x Red
ACC340B	1 x Blue
ACC340B-50	50 x Blue

Body Aid Ear Tips

Available in 3 sizes.
Supplied individually
or in packs of 50.



ACC263A	1 x Small
ACC263A-50	50 x Small
ACC263B	1 x Medium
ACC263B-50	50 x Medium
ACC263C	1 x Large
ACC263C-50	50 x Large

Silicone Ear Tips

Available in 5 sizes.
Supplied in pairs.



ACC102L-BLACK	12mm - Short
ACC102L-BLUE	10.5mm - Long
ACC102L-GREEN	10.5mm - Short
ACC102L-ORANGE	14mm - Long
ACC102L-RED	12mm - Long

Foam Ear Tips

Available in 5 sizes.
Supplied in pairs.




ACC102M-BLACK	12mm - Short
ACC102M-BLUE	10.5mm - Long
ACC102M-GREEN	10.5mm - Short
ACC102M-ORANGE	14mm - Long
ACC102M-RED	12mm - Long

Pre-Assembled Tubing Assemblies

Featuring eartip, tubing and adaptor. Available in 5 sizes.



Invisi Loop
Straight



Invisi Loop
Coiled



Coiled Tubing
Assembly

ACC396A	Invisi Loop Straight - Size 1
ACC396B	Invisi Loop Straight - Size 2
ACC396C	Invisi Loop Straight - Size 3
ACC396D	Invisi Loop Straight - Size 4
ACC396E	Invisi Loop Straight - Size 5

ACC397A	Invisi Loop Coiled - Size 1
ACC397B	Invisi Loop Coiled - Size 2
ACC397C	Invisi Loop Coiled - Size 3
ACC397D	Invisi Loop Coiled - Size 4
ACC397E	Invisi Loop Coiled - Size 5

ACC148N	Coiled Tubing Assembly
----------------	------------------------

Tubing, Tips, Earmoulds & Hooks

Generic Fit Earmoulds

Soft pliable open earmould attenuator. Available in 3 sizes for left or right ears, with or without link tubing. Presenters' version also available with coiled tubing, attenuator and standard earphone adaptor. Presenters' earmould replacement tubing also available.

Standard Earmould - Without Tubing

ACC149S	Small - Left
ACC149T	Medium - Left
ACC149U	Large - Left
ACC149V	Small - Right
ACC149W	Medium - Right
ACC149X	Large - Right

Standard Earmould - With Tubing

ACC149L	Small - Left
ACC149M	Medium - Left
ACC149N	Large - Left
ACC149P	Small - Right
ACC149Q	Medium - Right
ACC149R	Large - Right

Standard Presenters' Earmould - With Coiled Tubing, Attenuator & Earphone Adaptor


ACC342A	Small - Left
ACC342B	Medium - Left
ACC342C	Large - Left
ACC342D	Small - Right
ACC342E	Medium - Right
ACC342F	Large - Right

Standard Presenters' Earmould - 90° Bend - With Coiled Tubing, Attenuator & Earphone Adaptor

ACC342H	Small - Left
ACC342J	Medium - Left
ACC342K	Large - Left
ACC342L	Small - Right
ACC342N	Medium - Right
ACC342M	Large - Right

Generic Presenters' Earmould

ACC148H-SML-L	Small - Left
ACC148H-MED-L	Medium - Left
ACC148H-LRG-L	Large - Left
ACC148H-SML-R	Small - Right
ACC148H-MED-R	Medium - Right
ACC148H-LRG-R	Large - Right



Shirt Clips For Presenters Moulds

ACC148L	Shirt Clip
ACC148LNC	Broadcast Quality Shirt Clip

Receiver Adaptor With Metal Clip

ACC145D	Adaptor With Metal Clip
---------	-------------------------

Earmould Attenuator

ACC342G	Earmould Attenuator
---------	---------------------

Presenters' Earmould Replacement Tubing

ACC148RTC-COM	Complete Coiled Tube
ACC148RTH	Hard Tube
ACC148RTS	Straight Tube
ACC148RTS-COM	Complete Straight Tube

Coiled Tubing

ACC148Q	Clear - Left Or Right Ear
ACC148Q-BLK	Black - Left Or Right Ear
ACC148Q-DISL	DisappEar - Left
ACC148Q-DISR	DisappEar - Right

Tubing, Tips, Earmoulds & Hooks

General Earhooks

Replacement hooks for behind-the-ear hearing systems.



Style 1



Style 2



Style 3



Style 4



Style 5

ACC103B Style 1 - Standard

ACC105A Style 1 - Filtered

ACC104 Style 2

ACC105B Style 3 - Filtered

ACC103A Style 4

ACC271 Style 5

Puretone BTE Earhooks

Replacement hooks for Puretone behind-the-ear hearing systems.



PT-4/PT-7/PT-70

ACC257H



PT-80

ACC257J



PT-90

ACC257K



Xenon XS/Krypton-32

DAC-HOOK-001



Neon/Xenon

INT020

ACC257H PT-4/PT-7/PT-70

ACC257J PT-80

ACC257K PT-90

DAC-HOOK-001 Xenon XS/Krypton-32

INT020 Neon/Xenon - Standard

INT021 Neon/Xenon - Filtered

Tubing, Tips, Earmoulds & Hooks

Puretone P Series BTE Earhooks



Replacement hooks for Puretone P Series behind-the-ear hearing systems.



PS-HOOK1 P640/P660SP/P660 - Standard

PS-HOOK2 P680MW - Standard

PS-HOOK3 P680SW/P1200SW/P1200MW
Standard

PS-HOOK4 P660SP - Hook With Filter

PS-MHOOK1 P640/P660SP/P660 - Mini Hook

PS-MHOOK2 P680MW - Mini Hook

PS-MHOOK3 P680SW/P1200SW/P1200MW
Mini Hook

Puretone BTE Open Tube Adaptor

For use with 22MA, PT-40 and PT-70.



PT-BTE-OTA Open Tube Adaptor

P Series Open Tube Adaptor



PS-OTA1 for use with P640S/P660M.

PS-OTA2 for use with P680SW/P680M.



PS-OTA1 Open Tube Adaptor - P640S/P660M

PS-OTA2 Open Tube Adaptor - P680SW/P680MW